



NETS

SAFELY REMOVE FOREIGN BODIES AND POLYPS

The MT net has been developed for the removal of polyps and foreign bodies and enables the reliable recovery of tissue samples. The net itself consists of a reinforced thread material, which offers additional safety during recovery. Thanks to its unique dimensional stability and large opening, it facilitates work over several passes.

Twister® Plus is a 360° rotatable net that allows polyps or foreign objects to be easily and safely removed. In the large, three-dimensional mesh, several polyps and foreign objects can be detected simultaneously in one pass. The material is tear-resistant and of high quality. Four net sizes and three system lengths are available. Twister® Plus thus offers the ideal solution for every requirement.

SPECIFIC CHARACTERISTICS

- Large net
- High quality manufacture
- Comfortable 3-ring handle
- Easy and precise opening



Twister® Plus in different net sizes

Shape stable MT net

SPECIFICATIONS

REF	Ø net mm	Ø catheter mm	Working length mm	Ø working channel mm
TWISTER® PLUS				
HIC-GDL-27-5	22	2.5	230	2.8
HIC-GDL-26-5	26	2.5	230	2.8
HIC-GDL-29-5	32	2.5	230	2.8
HIC-GDL-30-5	15	1.7	120	2.0
HIC-GDL-31-5	15	1.7	160	2.0

Packaging unit: 5 pieces

REF	Ø net mm	Snare	Working length mm	Ø working channel mm
RETRIEVAL NET				
NT-E-30/23-P	30	Oval	2300	2.8
NT-E-30/16-P	30	Oval	1600	2.8

Packaging unit: 2 pieces