



MICRO-TECHTM
ENDOSCOPY

SINGLE-USE HEMOCLIP FAMILY

SURECLIPTM

Repositioning
Reliable 360° Rotation
Shorter Stem



YOUR TRUSTED PARTNER IN ENDOSCOPY

CLIP WITH CONFIDENCE

REPOSITIONING

SureClip's unique design permits repositioning prior to deployment. Being able to reposition a clip, may help improve placement accuracy.

RELIABLE 360° ROTATION

SureClip can be rotated, which helps to provide the correct orientation for tissue approximation or vessel closure.

SHORTER STEM

A shorter stem makes the clip less obtrusive, improving visualization of the target area, particularly when multiple clips are placed in close proximity.



SPECIFICATIONS

	REF	Clip Opening(mm)	Working Length(cm)	Min Working Channel(mm)
SureClip™ mini	ROCC-C-26-165-C	8	165	2.8
	ROCC-C-26-195-C	8	195	2.8
	ROCC-C-26-230-C	8	230	2.8
SureClip™	ROCC-D-26-165-C	11	165	2.8
	ROCC-D-26-195-C	11	195	2.8
	ROCC-D-26-230-C	11	230	2.8
SureClip™ plus	ROCC-F-26-165-C	16	165	2.8
	ROCC-F-26-195-C	16	195	2.8
	ROCC-F-26-230-C	16	230	2.8

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.

YOUR TRUSTED PARTNER IN ENDOSCOPY

MICRO-TECH (U.K.) Ltd | The Guild Hall, St Martins Courtyard, Coney Street, York, YO1 9QL
P: +44 (0)20 8016 1990 | sales@micro-tech-uk.com | www.micro-tech-uk.com